



COMPLEX HEALTHY SKIN

WHY IS SKIN HEALTH IMPORTANT?



The skin is the outer covering of the body, actively participating in all its life processes, performing a number of important functions: respiratory, protective, secretory, receptor, and others. It is the largest organ of the human body, often one of the first to signal internal organ diseases.



With age, a person's adaptive capabilities of all organs and tissues, including the skin, deteriorate. In elderly and old age, people's skin is more prone to various dermatoses, dryness, and the appearance of spots. During this period, any damage to the skin can have a long period of regeneration and often becomes infected.



Skin damage largely depends on the phototype - the lighter the skin, the more susceptible it is to the negative effects of environmental factors (sun, cold, wind), and the faster it loses elasticity, ceases to perform protective functions, becomes thin, and develops pigmentation.



Failure to observe personal hygiene or weakening of protective functions can allow pathogenic fungi to enter the body. As a rule, this occurs through personal contact with an infected person, their belongings, or by walking on surfaces contaminated with fungi (baths, floors). Pathogens are capable of affecting not only the skin, nails, and hair but also underlying tissues and even organs.



DeVita AP – the action is based on the phenomenon of resonance, in particular, the studies of the scientist and physician Royal Raymond Rife. Rife discovered that each microorganism has an individual vibration frequency, an increase in the amplitude of which leads to the death of the microorganism without harming other cells.

Viral, bacterial, and fungal infections, harmful habits, poor nutrition and care, frequent stress are causes of skin, nail, and hair problems. The DeVita AP programs contribute to reducing inflammatory reactions, decreasing infectious and toxic load.

The “Cleansing” automatic program is used to cleanse the body after applying a complex of antiparasitic programs.



DeVita Ritm – here the action is based on the studies of the scientists Dr. Reinhold Voll and A.F. Cramer, which found that each organ has a healthy frequency of oscillation. DeVita Ritm gently affects the body, pointing out the frequency of the vibrations of healthy organs. As a result, the body is capable of self-healing.

The automatic programs “Healthy Day - automatic” and “Health - automatic” have a beneficial effect on the work of the whole organism, increasing its resistance to foreign agents. They strengthen the defenses, increase vitality. The “Complete Cleansing” automatic program promotes the restoration and harmonization of the body’s functioning.

Using special programs activates metabolic processes, improves blood microcirculation, alleviates inflammatory manifestations, improves the condition of the skin, nails, and hair, and slows down the aging process.



DeVita Energy harmonizes the body’s work, energizes, activates the removal of toxins, activates the processes of self-regulation.

DIRECTIONS FOR USE



NUTRITIONAL RECOMMENDATIONS FOR USING DEVITA

Do not drink carbonated alcoholic and non-alcoholic drinks - they reduce the effectiveness of the devices. Remove or minimize coffee in your diet.

Avoid eating “junk food” (chips, crackers, etc.) - this can cause digestive problems. Drink clean drinking water at least 1-1.5 L / day to remove toxins from the body.



USE DEVITA DEVICES IN THE FOLLOWING SITUATIONS:

- for inflammatory conditions on the skin and nails;
- occurrence of pigmentation, warts, cracks;
- skin damage.



COMBINE THE DEVICES WITH DELIXIR SMART CELL NUTRITION

For best results, it is recommended to use the DeVita devices in conjunction with the DeLixir Collagen+ and DeLixir Detox smart cell nutrition supplements.

SCHEME NO.1 FIRST MONTH	According to the scheme of 1st month of the «Clean skin» complex
----------------------------	------------------------------------------------------------------

SCHEME NO.2 SECOND MONTH	DEVICE	PROGRAM		APPLICATION	
	1	DeVita Energy	Basic	Energy or Energy Pro	Once in the morning
	2	DeVita Ritm	Automatic programs	«Healthy day»*, «Health»*	It the daytime after DeVita Energy programs or from 8:00. Alternating every other day
	3		Special programs	«Healthy skin» or «Smoothing of wrinkles» and «Facial tone» to maintain youthfulness	After DeVita Ritm «Healthy day» program
	4	DeVita AP	Automatic programs	One of the 3 automatic programs: «Skin without pimples», «Healthy nails», «Stop skin fungus»*	In the nightttime. Alternate one of 3 automatic programs along with Cleansing program every other day.
	5		Special programs	Any of the manual programs depending on the problem. «Without toxins»	If necessary after DeVita Ritm «Health» program.

SCHEME NO.3 THIRD MONTH	DEVICE	PROGRAM		APPLICATION	
	1	DeVita Ritm	Automatic programs	«Healthy day»*, «Complete cleansing»*	In the daytime from 8:00. Alternating every other day
	2		Special programs	«Healthy skin» or «Smoothing of wrinkles» and «Facial tone» to maintain youthfulness	After DeVita Ritm «Healthy day» program
3	DeVita AP	Automatic programs	Stop parasites*	In the nightttime daily for 1 week.	

SCHEME NO.4 FOURTH MONTH	DEVICE	PROGRAM		APPLICATION	
	1	DeVita Ritm	Automatic programs	«Healthy day»*, «Health»*	In the daytime from 8:00. Alternating every other day.
	2		Special programs	«Healthy skin» or «Smoothing of wrinkles» and «Facial tone» to maintain youthfulness	After DeVita Ritm «Healthy day» program
	3	DeVita AP	Automatic programs	One of the 3 automatic programs: «Skin without pimples», «Healthy nails», «Stop skin fungus»	At the nightttime for 1 month.
4	Special programs		Any of the manual programs depending on the problem. «Without toxins»	If necessary after DeVita Ritm «Health» program	

*Programs of «Clean skin» complex.

RECOMMENDATIONS TO MAINTAIN SKIN HEALTH



Promptly treat acute and chronic illnesses



Eat a rational and balanced diet.



Consume an adequate amount of fluids.



Eliminate unhealthy habits.



Follow a rational work and rest schedule, get enough sleep.



Avoid aggressive exposure of the environment to the skin.



Maintain personal hygiene: take a shower daily, do not use others' hygiene products and cosmetics, use antibacterial agents less frequently.



You will get the most effective results by combining the preventive measures with the use of DeVita technologies and DeLixir smart cell nutrition supplements

Before use



Inflammation, redness, scaling, itching.

After 2 days of using DeVita



Reduction of inflammation, itching.

After a month of using DeVita



Improvement of the skin, hair, and nails condition. Good mood and well-being.



Use the scheme of the "HEALTHY SKIN" complex with the DeLixir COLLAGEN+ and DeLixir DETOX smart cell nutrition supplements, and your results will be noticeable much faster