



COMPLEX MOVEMENT WITHOUT PAIN

WHY SPINE AND JOINT PROBLEMS ARE DANGEROUS?



Problems with spine and joints are a fee that's humans have to pay for upright posture. Having only two points of support, the human body is forced to constantly maintain balance. Strictly speaking, each person is fated with problems with the spine. Because of this, the joints experience increased stress.



We exacerbate these problems by leading a sedentary lifestyle, lifting weights, or not paying attention to posture. A huge role is played by shoes, which are designed to soften the load on the spine, but do not always fulfill this function. Meanwhile, over 70% of people have spinal curvature.



The most common problems are impaired mobility, severe pain and inflammation. And this indirectly affects the appearance, causing obesity, varicose veins, problems with gastrointestinal tract and general fatigue. In especially severe cases, problems with the spine and joints lead to restricted movement and disability.



Diseases of the back and joints can go unnoticed for a long time, but at the most unexpected moment declare about themselves with severe pain. Treatment can, basically, only relieve symptoms - reduce pain, improve mobility and blood circulation. Therefore, it is important to take care of the health of the spine and joints on time.



DeVita AP – the action is based on the phenomenon of resonance, in particular, the studies of the scientist and physician Royal Raymond Rife. Rife discovered that each microorganism has an individual vibration frequency, an increase in the amplitude of which leads to the death of the microorganism without harming other cells.

The parasite larvae that enter the body settle on the membrane lining the joint cavity, or in the periarticular tissues, causing inflammation. Parasites also affect the spine. DeVita AP basic antiparasitic programs affect most known pathogens, quickly, effectively and without side effects, removing the parasitic load from the organism.

“Cleansing” automatic program is used to cleanse the body after applying basic antiparasitic programs, reduce the parasitic intoxication and improve overall well-being.

Special programs “No bacteria” and “Drainage Antitox” help to reduce toxic load in the presence of specific groups of microorganisms and salts of heavy metals in the body that have a negative effect on the musculoskeletal system.



DeVita Ritm – In this device, the development of scientists R. Voll and A. F. Kramer is embodied, which found that each organ has a healthy frequency of oscillations. DeVita Ritm gently affects the body, by giving it the frequency of vibrations of healthy organs. As a result, the body works on self-healing.

“Healthy Day” and «Health» automatic programs improve blood circulation, gently reduce the toxic load, which enhances immunity, improve the outflow of metabolic products from damaged areas. “Complete Cleansing” automatic program helps in main excretory system restoration and harmonizes the organism.

The programs “Recovery from injuries”, “Healthy joints - automatic” and “Joints without pain - automatic” relieve from inflammation, pain, swelling, activate metabolic processes in damaged tissues, normalize muscle tone, improve the functioning of the ligamentous apparatus.



DeVita Energy triggers self-regulation mechanisms, restores physiological processes in organs and tissues. Gives vigor and energy for movement.

DIRECTIONS FOR USE



NUTRITIONAL RECOMMENDATIONS FOR USING DEVITA

Do not drink carbonated alcoholic and non-alcoholic drinks - they reduce the impact of the devices. Remove or minimize coffee in your diet. Avoid eating “junk food” (chips, crackers, etc.) - this can cause digestive problems. Drink clean drinking water, at least 1-1.5 L / day to remove toxins from the body.



USE DEVITA DEVICES IN THE FOLLOWING SITUATIONS:

- pain in the joints and spine;
- injuries and inflammation of the muscles and joints;
- pinched nerves.



COMBINE THE DEVICES WITH DELIXIR SMART CELL NUTRITION

For best results, it is recommended to use the DeVita devices in conjunction with the DeLixir MULTI-ENERGY, DeLixir COLLAGEN+ and DeLixir DETOX smart cell nutrition supplements.









SCHEME №1 DURING 1 st MONTH	DEVICE	PROGRAM		USAGE	
	1	DeVita Energy	Energy or Energy Pro		Once in the morning
	2	DeVita Ritm	Auto-programs	Healthy Day-automatic, Health-automatic	Daily in the morning, alternating with each other every other day.
	3		Special programs	Joints without pain - automatic or Healthy joints and / or Recovery from injuries	Daily in the afternoon, if necessary.
	4	DeVita AP	Auto-programs	Mode 1 automatic	Daily at night during the 1 st week.
	5			Mode 2 automatic	Daily at night during the 2 nd week.
	6			Mode 3 automatic	Daily at night during the 3 rd week.
	7			Mode 4 automatic	Daily at night during the 4 th week.

SCHEME №2 DURING 2 nd MONTH	DEVICE	PROGRAM		USAGE	
	1	DeVita Energy	Energy or Energy Pro		Once in the morning
	2	DeVita Ritm	Auto-programs	Healthy Day-automatic, Health-automatic	Daily in the morning, alternating with each other every other day.
	3		Special programs	Joints without pain - automatic or Healthy joints and / or Recovery from injuries	Daily in the afternoon, if necessary.
	4	DeVita AP	Auto-programs	Cleansing-automatic	Daily at night
5	Special programs		No bacteria, Drainage Antitox	Daily in the evening after DeVita Ritm programs for 14 days.	

SCHEME №3 DURING 3 rd MONTH	DEVICE	PROGRAM		USAGE	
	1	DeVita Energy	Energy or Energy Pro		Once in the morning
	2	DeVita Ritm	Auto-programs	Complete Cleansing-automatic	Daily from 8 a.m.
3	DeVita AP	Auto-programs	Mode 1-4 automatic	If needed, at night.	

SCHEME №4 DURING 4 th MONTH	DEVICE	PROGRAM		USAGE	
	1	DeVita Energy	Energy or Energy Pro		Once in the morning
	2	DeVita Ritm	Auto-programs	Healthy Day-automatic, Health-automatic	Daily in the morning, alternating with each other every other day.
	3		Special programs	Joints without pain - automatic or Healthy joints and / or Recovery from injuries	Daily in the afternoon, if necessary.
4	DeVita AP	Special programs	No nematodes and ascarides, No parasites - automatic	At night, alternating every other day for 2 weeks	

PREVENTIVE MEASURES FOR SPINAL AND JOINT PROBLEMS

-  Use vitamins and minerals in spring and autumn.
-  Watch your posture. Lean your back against a flat wall - your nape, shoulder blades, buttocks and feet should touch the wall. If necessary, consult a doctor and purchase a corrective corset.
-  Sleep on an orthopedic mattress and pillow. In extreme cases - on a dense hard surface. Fall asleep lying on your back.
-  Exercise and strengthen your back muscles and abs. Developed abdominal muscles relieve pressure on the lower back.
-  When doing sports and lifting weights, follow the correct technique - do not round your back, strain your abs.
-  Eat sour-milk products rich in calcium, as well as products with gelatin content – meat with gelatin and jelly.
-  Consult a doctor if you want to try a new physical activity - yoga, step aerobics, strength training, etc.
-  Be wary of osteopathy, manual therapy, mechanical massage.



You will get the most effective results by combining the preventive measures with the use of DeVita technologies and DeLixir smart cell nutrition supplements.

Before
use



Pain in the back and joints,
fatigue, headache.

After the first use
of DeVita



Decreased pain
and tiredness.

After a month
of using DeVita



Fast recovery,
ease of motion,
no pain.



Use the scheme of the "MOVEMENT WITHOUT PAIN" complex together with the DeLixir MULTI-ENERGY, DeLixir COLLAGEN+ and DeLixir DETOX smart cell nutrition supplements, and your results will be noticeable much faster.